

Instruction Guide in Activating Your Megamatcher Dongle (Volume License)

Thank you for your purchase of the Megamatcher Dongle License. With this license, your company can use the Megamatcher Software, which allows fast fingerprint matching using the Department of Health IDTOMIS software for your drug testing operations.



<http://www.chiralssoft.com/megamatcher.php>

Use and Handling

1. Make sure your computer is connected to a reliable Internet connection.
2. Plug in the USB Dongle to a USB Computer Port preferably at the Back port of the computer to avoid interruption while the machine is being used.
3. Open Activation Wizard


Path: C:\Program Files\Neurotechnologija\Activation\ActivationWizardDotNet.exe

OR Start->Programs->Neurotechnologija->Activation Wizard


4. Select Volume License, click Next




5. In the first box, select this computer is “USING DONGLE”



6. Click Install License Manager



A window will appear like this



7. To check if the installation was successful, you can click on the “Diagnostic” button and look for the indication that “IsReg” is in blue color font.
8. You may also test by using IDITOMIS and taking fingerprints. Make sure your biometric fingerprint scanner is plugged into a USB port on your computer. If you are able to take fingerprints, your Megamatcher Dongle license is in working order.
9. Please note, every time you need to take fingerprints, your Dongle has to be plugged into the USB port the whole time, otherwise your biometric fingerprint scanner will not be able to take fingerprints.